

# **Growth and Maintenance of Flp-In<sup>™</sup> Cell Lines**

**Catalog nos. R750-07, R752-07, R758-07, R760-07, R761-07, R762-07**

**Version E**

February 12, 2003

25-0323



**[www.invitrogen.com](http://www.invitrogen.com)**  
**[tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)**



# Table of Contents

Table of Contents .....	iii
Important Information .....	iv
Accessory Products.....	v
Overview .....	1
<b>Methods .....</b>	<b>5</b>
Culturing Flp-In <sup>™</sup> Cell Lines .....	5
Freezing Cells .....	8
Transfecting Cells .....	9
<b>Appendix.....</b>	<b>12</b>
Technical Service.....	12
Purchaser Notification .....	14
References .....	16

# Important Information

---

## Shipping/Storage

All cell lines are shipped on dry ice. Store in liquid nitrogen upon receipt.

---

## Contents

This manual is supplied with the following cell lines:

Cell Line	Catalog no.
Flp-In <sup>™</sup> -293	R750-07
Flp-In <sup>™</sup> -CV-1	R752-07
Flp-In <sup>™</sup> -CHO	R758-07
Flp-In <sup>™</sup> -BHK	R760-07
Flp-In <sup>™</sup> -3T3	R761-07
Flp-In <sup>™</sup> -Jurkat	R762-07

All cell lines are supplied as one vial containing  $3 \times 10^6$  frozen cells in 1 ml of Freezing Medium.

---

## Product Qualification

The following criteria are used to qualify the Flp-In<sup>™</sup> cell lines:

- Each cell line is tested independently and certified to be free of mycoplasma.
  - Prior to freezing, cells are greater than 95% viable. Forty-eight hours after thawing, cells are greater than 90% viable.
  - Each cell line is tested for  $\beta$ -galactosidase activity by plating cells into 6-well plates in medium containing Zeocin<sup>™</sup>. After several days of growth, the cells are assayed for  $\beta$ -galactosidase activity using the  $\beta$ -Gal Staining Kit from Invitrogen. Each cell line must exhibit greater than 95% cells expressing  $\beta$ -galactosidase.
-

## Accessory Products

### Introduction

The products listed in this section are intended for use with the Flp-In™ Cell Lines and the Flp-In™ System. For more information, refer to our Web site ([www.invitrogen.com](http://www.invitrogen.com)) or call Technical Service (see page 12).

### Cell Culture Reagents

A large variety of Gibco™ cell culture products are available from Invitrogen to facilitate growth and maintenance of the Flp-In™ cell lines. For more information about the products listed below, refer to our Web site ([www.invitrogen.com](http://www.invitrogen.com)) or call Technical Service (see page 12). **Note:** Reagents are available in other sizes.

Item	Amount	Catalog no.
Dulbecco's Modified Eagle Medium (D-MEM)	500 ml	11965-092
Ham's F-12	500 ml	11765-054
RPMI Medium 1640	500 ml	11875-093
Fetal Bovine Serum	500 ml	16000-044
Donor Calf Serum	500 ml	16030-074
200 mM L-Glutamine	100 ml	25030-081
Penicillin-Streptomycin	100 ml	15070-063
Trypsin-EDTA	100 ml	25300-054

### Additional Reagents

The products listed below may be used with the Flp-In™ Cell Lines. Zeocin™ is available for maintenance and growth of the Flp-In™ cell lines. For more information, refer to our Web site ([www.invitrogen.com](http://www.invitrogen.com)) or call Technical Service (see page 12).

Item	Amount	Catalog no.
Zeocin™	1 g	R250-01
	5 g	R250-05
Hygromycin	1 g	R220-05
Lipofectamine™ 2000 Reagent	0.75 ml	11668-027
	1.5 ml	11668-019
Lipofectamine™ Reagent	1 ml	18324-012
Plus™ Reagent	0.85 ml	11514-015
Phosphate-Buffered Saline (PBS) 7.4 (1X)	500 ml	10010-023

*continued on next page*

## Accessory Products, continued

### Flp-In™ Products

The plasmids required to generate Flp-In™ host cell lines and expression cell lines are available separately from Invitrogen. For more information about the features of each vector, refer to our Web site ([www.invitrogen.com](http://www.invitrogen.com)) or call Technical Service (see page 12). Ordering information is provided below.

Product	Amount	Catalog no.
pFRT/ <i>lacZeo</i>	20 µg, lyophilized in TE	V6015-20
pFRT/ <i>lacZeo2</i>	20 µg, lyophilized in TE	V6022-20
pOG44	20 µg, lyophilized in TE	V6005-20
pcDNA™5/FRT	20 µg, lyophilized in TE	V6010-20
pcDNA™5/FRT/V5-His TOPO® TA Expression Kit	20 reactions	K6020-01
pSecTag/FRT/V5-His TOPO® TA Expression Kit	20 reactions	K6025-01
pEF5/FRT/V5-DEST Gateway® Vector Pack	6 µg	V6020-20

# Overview

## Introduction

The Flp-In™ cell lines stably express the *lacZ*-Zeocin™ fusion gene and are designed for use with the Flp-In™ System (Catalog nos. K6010-01 and K6010-02). Each cell line contains a single integrated Flp Recombination Target (FRT) site from pFRT/*lacZeo* or pFRT/*lacZeo2* as confirmed by Southern blot analysis. See below and the next page for information about the generation of the Flp-In™ cell lines. For more information about the Flp-In™ System and its components, refer to the Flp-In™ System manual, visit our Web site ([www.invitrogen.com](http://www.invitrogen.com)), or call Technical Service (see page 12). The Flp-In™ System manual is also available for downloading from our Web site.

Generation of Flp-In™ expression cell lines requires cotransfection of the Flp-In™ cell line with a Flp-In™ expression vector containing your gene of interest and the Flp recombinase expression plasmid, pOG44 (O'Gorman *et al.*, 1991). Flp recombinase mediates insertion of your Flp-In™ expression construct into the genome at the integrated FRT site through site-specific DNA recombination (O'Gorman *et al.*, 1991; Sauer, 1994). Stable cell lines expressing your gene of interest from the Flp-In™ expression vector can be generated by selection using hygromycin B. For more information about FRT sites and Flp recombinase-mediated DNA recombination, refer to the Flp-In™ System manual.

## Parental Cell Lines

The table below provides a brief description of the source of the parental cell line used to generate each Flp-In™ cell line. The parental cell lines were obtained from the American Type Culture Collection (ATCC). The ATCC number for each cell line is included. For further information about the parental cell lines, refer to the ATCC Web site ([www.atcc.org](http://www.atcc.org)).

Cell Line	Characteristic	Source	ATCC Number
293	Adherent	Human embryonic kidney (Graham <i>et al.</i> , 1977)	CRL-1573
CV-1	Adherent	African Green Monkey kidney (Kit <i>et al.</i> , 1965)	CCL-70
CHO-K1	Adherent	Chinese Hamster ovary (Kao and Puck, 1968)	CCL-61
BHK	Adherent	Baby hamster kidney (Talavera and Basilico, 1977)	CCL-10
NIH/3T3	Adherent	Mouse (NIH Swiss) embryonic fibroblast (Jainchill <i>et al.</i> , 1969)	CRL-1658
Jurkat	Suspension	Human T-cell leukemia (Weiss <i>et al.</i> , 1984)	TIB-152


*continued on next page*

## Overview, continued

### Description of Flp-In™ Cell Lines


All of the Flp-In™ cell lines (except Flp-In™-CHO; see below) contain a single integrated FRT site and stably express the *lacZ*-Zeocin™ fusion gene from the pFRT/*lacZeo* plasmid under the control of the SV40 early promoter (see diagram below). The location of the FRT site in each Flp-In™ cell line has not been mapped, but is presumed to have integrated into a transcriptionally active genomic locus as determined by generation of a Flp-In™ expression cell line containing the pcDNA™5/FRT/CAT or pEF5/FRT/GW-CAT control plasmid. The Flp-In™ cell lines should be maintained in medium containing Zeocin™ (see the next page).

For more information about pFRT/*lacZeo*, refer to the Flp-In™ System manual. For information about pcDNA™5/FRT/CAT or pEF5/FRT/GW-CAT, refer to the pcDNA™5/FRT or pEF5/FRT/V5-DEST manuals, respectively.



### Flp-In™-CHO Cell Line

The Flp-In™-CHO cell line contains a single integrated FRT site and stably expresses the *lacZ*-Zeocin™ fusion gene from the pFRT/*lacZeo*2 plasmid. Note that pFRT/*lacZeo*2 contains a mutated SV40 early promoter (P<sub>SV40Δ</sub>) which is severely abrogated in its activity. The SV40Δ early promoter in pFRT/*lacZeo*2 exhibits approximately 60-fold less activity than the wild-type SV40 early promoter in pFRT/*lacZeo*. Because of the minimal activity of the SV40Δ promoter, we expect that stable transfectants expressing the *lacZ*-Zeocin™ gene from pFRT/*lacZeo*2 should contain FRT sites which have integrated into the most transcriptionally active genomic loci. The location of the FRT site in the Flp-In™-CHO cell line has not been mapped, but has been demonstrated to have integrated into a highly transcriptionally active genomic locus as determined by generation of a Flp-In™ expression cell line containing the pcDNA™5/FRT/luc (luciferase-expressing) control plasmid. The Flp-In™-CHO cell line should be maintained in medium containing Zeocin™ (see below). For more information about pFRT/*lacZeo*2 and the SV40Δ early promoter, refer to the pFRT/*lacZeo*2 manual.



*continued on next page*



## Overview, continued

### Media for Cell Lines

The table below lists the recommended complete medium, freezing medium, and antibiotic concentration required to maintain and culture each Flp-In™ cell line.

Cell Line	Complete Medium	[Antibiotic]	Freezing Medium
Flp-In™-293	D-MEM (high glucose) 10% FBS* 2 mM L-glutamine 1% Pen-Strep (optional)	100 µg/ml Zeocin™	90% complete medium 10% DMSO
Flp-In™-CV-1	D-MEM (high glucose) 10% FBS* 2 mM L-glutamine 1% Pen-Strep (optional)	100 µg/ml Zeocin™	90% complete medium 10% DMSO
Flp-In™-CHO	Ham's F12 10% FBS* 2 mM L-glutamine 1% Pen-Strep (optional)	100 µg/ml Zeocin™	90% complete medium 10% DMSO
Flp-In™-BHK	D-MEM (high glucose) 10% FBS* 2 mM L-glutamine 1% Pen-Strep (optional)	100 µg/ml Zeocin™	90% complete medium 10% DMSO
Flp-In™-3T3	D-MEM (high glucose) 10% donor calf serum 2 mM L-glutamine 1% Pen-Strep (optional)	100 µg/ml Zeocin™	90% complete medium 10% DMSO
Flp-In™-Jurkat	RPMI 1640 10% FBS* 2 mM L-glutamine 1% Pen-Strep (optional)	100 µg/ml Zeocin™	90% complete medium 10% DMSO

\*FBS = fetal bovine serum

*continued on next page*

## Overview, continued

### Important Guidelines

- 
- FBS does not need to be heat inactivated for use with these cell lines.
  - Cell lines should be maintained in medium containing Zeocin™ at the concentrations listed above.
  - If adherent cells (e.g. Flp-In™-293, Flp-In™-CV-1, Flp-In™-CHO, Flp-In™-3T3, Flp-In™-BHK) are split at a 1:5 to 1:10 dilution, they will generally reach 80-90% confluence in 3-4 days.
  - Suspension Flp-In™-Jurkat cells will demonstrate optimal growth characteristics if maintained at a cell density between  $1 \times 10^5$  cells/ml and  $1 \times 10^6$  cells/ml.
  - When maintaining Flp-In™-Jurkat cells in suspension culture, do not allow the medium to turn yellow; this indicates that cells have reached too high a density or that the medium is depleted of nutrients. If this occurs, either add fresh complete media to the cells or passage them.
-

# Methods

## Culturing Flp-In™ Cell Lines

---

### General Cell Handling

Follow the guidelines below to successfully grow and maintain your cells.

- All solutions and equipment that come in contact with the cells must be sterile. Always use proper sterile technique and work in a laminar flow hood.
- Before starting experiments, be sure to have cells established and also have some frozen stocks on hand. We recommend that you always use early-passage cells for your experiments. Upon receipt of the cells from Invitrogen, grow and freeze multiple vials of the particular cell line to ensure that you have any adequate supply of early-passage cells.
- Cells should be at the appropriate confluence (approximately 60%) and >90% viability prior to transfection (see page 9).
- For general maintenance of cells, pass all cell lines when they are 80-90% confluent (for adherent cells) or when they reach a density of  $2 \times 10^6$  to  $4 \times 10^6$  cells/ml (for suspension cells).
- Use trypan blue exclusion to determine cell viability. Log phase cultures should be >90% viable.

---

### Before Starting

Be sure to have the following solutions and supplies available:

- 15 ml sterile, conical tubes
- 5, 10, and 25 ml sterile pipettes
- Cryovials
- Phosphate-Buffered Saline (see page v for ordering information)
- 0.4% Trypan blue in PBS and hemacytometer (for counting cells)
- Reagents to prepare the appropriate complete medium (see page 3)
- Freezing Medium (see pages 3 and 8)
- Table-top centrifuge
- 75 cm<sup>2</sup> flasks, 175 cm<sup>2</sup> flasks and other appropriately-sized tissue culture flasks or plates
- Trypsin/versene (EDTA) solution or other trypsin solution

---

*continued on next page*

## Culturing Flp-In™ Cell Lines, continued

---

### Thawing Adherent Cells

The following protocol is designed to help you thaw adherent cells to initiate cell culture. All cell lines are supplied in vials containing  $3 \times 10^6$  cells in 1 ml of Freezing Medium.

1. Remove the vial of cells from the liquid nitrogen and thaw quickly in a 37°C water bath.
  2. Just before the cells are completely thawed, decontaminate the outside of the vial with 70% ethanol, and transfer the cells to a T-75 flask containing 12 ml of complete medium without Zeocin™.
  3. Incubate the flask at 37°C for 2-4 hours to allow the cells to attach to the bottom of the flask.
  4. Aspirate off the medium and replace with 12 ml of fresh, complete medium without Zeocin™.
  5. Incubate cells overnight at 37°C.
  6. The next day, aspirate off the medium and replace with fresh, complete medium containing Zeocin™ (at the recommended concentration listed on page 3).
  7. Incubate the cells and check them daily until the cells are 80-90% confluent (2-7 days).
  8. Proceed to **Passaging Adherent Cells**, next page.
- 

### Thawing Suspension Cells

The following protocol is designed to help you thaw suspension cells to initiate cell culture. All cell lines are supplied in vials containing  $3 \times 10^6$  cells in 1 ml of Freezing Medium.

1. Remove the vial of cells from the liquid nitrogen and thaw quickly in a 37°C water bath.
  2. Just before the cells are completely thawed, decontaminate the outside of the vial with 70% ethanol, and transfer the cells to a sterile, conical tube containing 5 ml of complete medium without Zeocin™.
  3. Centrifuge for 3 minutes at 750 x g at room temperature.
  4. Aspirate off the medium and resuspend the cells in 12 ml of fresh, complete medium without Zeocin™.
  5. Transfer the cells to a T-75 flask and incubate cells overnight at 37°C.
  6. The next day, add Zeocin™ to the cells (at the recommended concentration listed on page 3).
  7. Incubate the cells and count them daily until the cells reach a density ranging from  $2 \times 10^6$  cells/ml to  $4 \times 10^6$  cells/ml (2-7 days). **Note:** You may add fresh, complete medium containing Zeocin™ to the cells every few days.
  8. Proceed to **Passaging Suspension Cells**, next page.
- 

*continued on next page*

## Culturing Flp-In<sup>™</sup> Cell Lines, continued

---

### Passaging Adherent Cells

1. When cells are ~80-90% confluent, remove all medium from the flask.
  2. Wash cells once with 10 ml PBS to remove excess medium and serum. Serum contains inhibitors of trypsin.
  3. Add 5 ml of trypsin/versene (EDTA) solution to the monolayer and incubate 1 to 5 minutes at room temperature until cells detach. Check the cells under a microscope and confirm that most of the cells have detached. If cells are still attached, incubate a little longer until most of the cells have detached.
  4. Add 5 ml of complete medium to stop trypsinization.
  5. Briefly pipet the solution up and down to break up clumps of cells.
  6. To maintain cells in 75 cm<sup>2</sup> flasks, transfer 1 ml of the 10 ml cell suspension from Step 5 to a new 75 cm<sup>2</sup> flask and add 15 ml fresh, complete containing Zeocin<sup>™</sup>. If you want the cells to reach confluency sooner, split the cells at a lower dilution (e.g. 1:4).  
**Note:** To expand cells into 175 cm<sup>2</sup> flasks, add 28 ml of fresh, complete medium containing Zeocin<sup>™</sup> to each of three 175 cm<sup>2</sup> flasks, then transfer 2 ml of the cell suspension to each flask to obtain a total volume of 30 ml.
  7. Incubate flasks in a humidified, 37°C, 5% CO<sub>2</sub> incubator.
  8. Repeat Steps 1-7 as necessary to maintain or expand cells.
- 

### Passaging Suspension Cells

1. Passage suspension cells when they reach a density of 2-4 x 10<sup>6</sup> cells/ml.
  2. To maintain cells in 75 cm<sup>2</sup> flasks, transfer 1-1.5 ml of cell suspension from Step 1 to a new 75 cm<sup>2</sup> flask containing 13-14 ml of fresh, complete medium with Zeocin<sup>™</sup>.  
**Note:** You may split the cells at a lower dilution (e.g. 1:4), if desired.
  3. To expand cells into 175 cm<sup>2</sup> flasks, add 28 ml of fresh, complete medium containing Zeocin<sup>™</sup> to each of three 175 cm<sup>2</sup> flasks, then transfer 2 ml of the cell suspension to each flask to obtain a total volume of 30 ml.  
**Note:** You may also expand cells into a spinner flask, if desired.
  4. Incubate flasks in a humidified, 37°C, CO<sub>2</sub> incubator.
- Repeat Steps 1-4 as necessary to maintain or expand cells.
-

# Freezing Cells

---

## Introduction

When freezing the Flp-In™ cell lines, we recommend the following:

- Freeze cells at a density of **at least**  $3 \times 10^6$  cells/ml.
- Use a freezing medium composed of 90% complete medium and 10% DMSO. Complete medium is medium containing serum.

Guidelines to prepare freezing medium and freeze cells are provided in this section.

---

## Preparing Freezing Medium

Freezing medium should be prepared fresh immediately before use.

1. In a sterile, conical centrifuge tube, mix together the following reagents for every 1 ml of freezing medium needed:

Fresh complete medium	0.9 ml
DMSO	0.1 ml
  2. Place the tube on ice. Discard any remaining freezing medium after use.
- 

## Freezing the Cells

Before starting, label cryovials and prepare freezing medium (see above). Keep the freezing medium on ice.

1. To collect cells, perform the following:
  - For adherent cells, follow Steps 1-5 of **Passaging Adherent Cells**, page 7.
  - For suspension cells, transfer cells to a sterile, conical centrifuge tube.
2. Count the cells.
3. Pellet cells at  $250 \times g$  for 5 minutes in a table top centrifuge at room temperature and carefully aspirate off the medium.
4. Resuspend the cells at a density of **at least**  $3 \times 10^6$  cells/ml in chilled freezing medium.
5. Place vials in a microcentrifuge rack and aliquot 1 ml of the cell suspension into each cryovial.
6. Freeze cells in an automated or manual, controlled-rate freezing apparatus following standard procedures. For ideal cryopreservation, the freezing rate should be a decrease of  $1^\circ\text{C}$  per minute.
7. Transfer vials to liquid nitrogen for long-term storage.

**Note:** You may check the viability and recovery of frozen cells 24 hours after storing cryovials in liquid nitrogen by following the procedure outlined in **Thawing Adherent Cells** or **Thawing Suspension Cells**, page 6, as appropriate.

---

# Transfecting Cells

## Introduction

To generate stable Flp-In™ expression cell lines, you will cotransfect your Flp-In™ expression construct and the pOG44 plasmid into the Flp-In™ cell line and select for stable transfectants using hygromycin. General guidelines and recommendations for transfection are provided in this section. We recommend that you read through this section before beginning.

## Transfection Methods

The Flp-In™ cell lines are generally amenable to transfection using standard methods including calcium phosphate precipitation (Chen and Okayama, 1987; Wigler *et al.*, 1977), lipid-mediated transfection (Felgner *et al.*, 1989; Felgner and Ringold, 1989), and electroporation (Chu *et al.*, 1987; Shigekawa and Dower, 1988). We typically use lipid-based transfection reagents to introduce Flp-In™ expression constructs into the Flp-In cell lines. The table below lists the recommended transfection reagent for each Flp-In™ cell line. Other transfection reagents may be suitable.

To transfect...	Use...
Flp-In™-293 Flp-In™-CV-1 Flp-In™-CHO Flp-In™-Jurkat	Lipofectamine™ 2000 Reagent
Flp-In™-3T3 Flp-In™-BHK	Lipofectamine™ Reagent and Plus™ Reagent

**Note:** Lipofectamine™ 2000 Reagent, Lipofectamine™ Reagent, and Plus™ Reagent are available from Invitrogen (see page v for ordering information).



### Important

We have observed down-regulation of the viral CMV promoter and subsequent loss of gene expression when pcDNA™5/FRT-based expression constructs are introduced into Flp-In™-3T3 or Flp-In™-BHK cells. This behavior is not observed with pEF5/FRT-based expression constructs. If you are generating Flp-In™ expression cell lines using the Flp-In™-3T3 or Flp-In™-BHK cell line, we recommend that you express your gene of interest from a pEF5/FRT-based plasmid (*e.g.* pEF5/FRT/V5-DEST).

*continued on next page*

## Transfecting Cells, continued

### Generating Stable Expression Cell Lines

To generate Flp-In™ expression cell lines, cotransfect your Flp-In™ expression construct and the pOG44 plasmid into the Flp-In™ cell line of choice, and select for stable transfectants using hygromycin B. Before transfection, you may want to test the sensitivity of the Flp-In™ cell line to hygromycin B to more accurately determine the hygromycin B concentration to use for selection. A suggested range of hygromycin B concentrations to use for selection of your Flp-In™ expression vector is listed below. For more information, refer to the Flp-In™ System manual. Hygromycin B may be obtained from Invitrogen (see page v for ordering information).

**Important:** Following cotransfection, your Flp-In™ expression clones should become sensitive to Zeocin™; therefore, your selection medium should NOT contain Zeocin™.

Cell Line (After Transfection with Flp-In™ Expression Vector)	Estimated Hygromycin B Concentration (µg/ml)
Flp-In™ -293	100-200
Flp-In™ -CV-1	100-200
Flp-In™ -CHO	500-600
Flp-In™ -BHK	100-200
Flp-In™ -3T3	100-200
Flp-In™ -Jurkat	200-400



#### Note

When transfecting Flp-In™ -CHO cells, we recommend following these guidelines:

- 48 hours after transfection, split the cells directly into medium containing the appropriate concentration of hygromycin.
- Split the cells such that they are no more than 25% confluent. If the cells are too dense, the antibiotic will not kill the cells. Antibiotics work best on actively dividing cells.

### Polyclonal Selection of Isogenic Cell Lines

Because every Flp-In™ cell line contains a single integrated FRT site, all of the hygromycin-resistant foci that you obtain after cotransfection with the Flp-In™ expression vector and pOG44 should be isogenic (*i.e.* the Flp-In™ expression vector should integrate into the same genomic locus in every clone; therefore, all clones should be identical). To obtain stable expression cell lines, you may perform “polyclonal” selection and screening of your hygromycin-resistant cells. After hygromycin selection, simply pool the hygromycin-resistant foci and screen the entire population of cells for the following phenotypes:

- Zeocin™ sensitivity
- Lack of β-galactosidase activity
- Expression of the gene of interest

*continued on next page*



## Transfecting Cells, continued

---

### Selection of Individual Cell Lines

If desired, single hygromycin-resistant foci can be isolated and expanded to generate individual clonal cell lines. To isolate individual clones, simply pick 5-20 hygromycin-resistant foci and expand the cells. You may verify that your FLP-In™ expression construct has integrated into the FRT site by testing each clone for Zeocin™ sensitivity and lack of  $\beta$ -galactosidase activity. Select those clones that are hygromycin-resistant, Zeocin™-sensitive, and lack  $\beta$ -galactosidase activity, and assay for expression of your gene of interest.

---



#### Note

Note that in rare instances, it is possible to generate a FLP-In™ expression cell line in which the FLP-In™ expression plasmid has undergone both FLP recombinase-mediated integration into the FRT site and random integration into a second genomic site. In this case, clones will still exhibit hygromycin resistance. To test for these “second site integrants”, transfect the cells with the pOG44 plasmid and select for Zeocin™ resistance. The FLP recombinase should mediate excision of the FLP-In™ expression plasmid at the FRT site and restore the *lacZ*-Zeocin™ fusion gene. The resulting cells should exhibit  $\beta$ -galactosidase activity, Zeocin™ resistance, and continued expression of the gene of interest. Alternatively, you may perform Southern blot analysis to identify second site integrants if suitable restriction enzymes are selected.

---

# Appendix

## Technical Service

### World Wide Web



Visit the Invitrogen Web Resource using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe® Acrobat® (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your Web browser (Internet Explorer 5.0 or newer or Netscape 4.0 or newer), then enter the following location (or URL):

**[www.invitrogen.com](http://www.invitrogen.com)**

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

### Contact Us

For more information or technical assistance, please call, write, fax, or email. Additional international offices are listed on our Web page ([www.invitrogen.com](http://www.invitrogen.com)).

#### Corporate Headquarters:

Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, CA 92008 USA  
Tel: 1 760 603 7200  
Tel (Toll Free): 1 800 955 6288  
Fax: 1 760 602 6500  
E-mail:  
[tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)

#### Japanese Headquarters:

Invitrogen Japan K.K.  
Nihonbashi Hama-Cho Park Bldg. 4F  
2-35-4, Hama-Cho, Nihonbashi  
Tel: 81 3 3663 7972  
Fax: 81 3 3663 8242  
E-mail: [jpinfo@invitrogen.com](mailto:jpinfo@invitrogen.com)

#### European Headquarters:

Invitrogen Ltd  
Inchinnan Business Park  
3 Fountain Drive  
Paisley PA4 9RF, UK  
Tel: +44 (0) 141 814 6100  
Tech Fax: +44 (0) 141 814 6117  
E-mail:  
[eurotech@invitrogen.com](mailto:eurotech@invitrogen.com)

### MSDS Requests

To request an MSDS, visit our Web site at [www.invitrogen.com](http://www.invitrogen.com). On the home page, go to 'Technical Resources', select 'MSDS', and follow instructions on the page.

*continued on next page*

## Technical Service, continued

---

### Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, please contact our Technical Service Representatives.

Invitrogen warrants that all of its products will perform according to the specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives.

**Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.**

---

# Purchaser Notification

---

## Introduction

Use of the Flp-In™ Cell Lines and the Flp-In™ System (“System”) is covered under a number of different licenses including those detailed below.

---

## Information for European Customers

The Flp-In™ cell lines are genetically modified and carry either the pUC-derived plasmid, pFRT/*lacZeo* or pFRT/*lacZeo2*. As a condition of sale, this product must be in accordance with all applicable local legislation and guidelines including EC Directive 90/219/EEC on the contained use of genetically modified organisms.

---

## Limited Use Label License No. 64: Flp-In™ System

Invitrogen Corporation (“Invitrogen”) has a license to sell the Flp-In™ System and its components (“System”) to scientists for research purposes only, under the terms described below. Use of the System for any Commercial Purpose (as defined below) requires the user to obtain commercial licenses as detailed below. Before using the System, please read the terms and conditions set forth below. Your use of the System shall constitute acknowledgment and acceptance of these terms and conditions. If you do not wish to use the System pursuant to these terms and conditions, please contact Invitrogen’s Technical Services within 10 days to return the unused and unopened System for a full refund. Otherwise, please complete the User Registration Card and return it to Invitrogen.

Invitrogen grants you a non-exclusive license to use the enclosed System for research purposes only. The System is being transferred to you in furtherance of, and reliance on, such license. You may not use the System, or the materials contained therein, for any Commercial Purpose without licenses for such purpose.

Commercial Purpose includes: any use of the System or Expression Products in a Commercial Product; any use of the System or Expression Products in the manufacture of a Commercial Product; any sale of the System or Expression Products; any use of the System or Expression Products to facilitate or advance research or development of a Commercial Product; and any use of the System or Expression Products to facilitate or advance any research or development program the results of which will be applied to the development of a Commercial Product.

“Expression Products” means products expressed with the System, or with the use of any vectors or host strains in the System. “Commercial Product” means any product intended for sale or commercial use.

Access to the System must be limited solely to those officers, employees and students of your entity who need access to perform the aforementioned research. Each such officer, employee and student must be informed of these terms and conditions and agree, in writing, to be bound by same. You may not distribute the System or the vectors or host strains contained in it to others. You may not transfer modified, altered, or original material from the System to a third party without written notification to, and written approval from Invitrogen. You may not assign, sub-license, rent, lease or otherwise transfer any of the rights or obligations set forth herein, except as expressly permitted by Invitrogen.

---

*continued on next page*

## Purchaser Notification, continued

---

**Limited Use Label  
License No. 64:  
Flp-In™ System,  
continued**

This product is licensed under U.S. Patent Nos. 5,654,182 and 5,677,177 and is for research purposes only. Inquiries about licensing for commercial or other uses should be directed to: The Salk Institute for Biological Studies, 10010 North Torrey Pines Road, La Jolla, CA 92037, Attn.: Department of Intellectual Property and Technology Transfer. Phone: 858-453-4100 ext 1275; Fax: 858-450-0509; Email: amueller@salk.edu.

---

**Limited Use Label  
License No. 86:  
Zeocin™ Selection  
Marker**

The *ble* resistance gene is the subject of one or more of U.S. Patent Nos. 5,021,344 and 5,118,620 and foreign equivalents owned by Cayla and licensed to Invitrogen. This product is sold for research purposes only. For commercial license information, please contact: Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, CA 92008. Phone: 760-603-7200, Fax: 760-602-6500.

---

## References

---

- Chen, C., and Okayama, H. (1987). High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. *Molec. Cell. Biol.* 7, 2745-2752.
- Chu, G., Hayakawa, H., and Berg, P. (1987). Electroporation for the Efficient Transfection of Mammalian Cells with DNA. *Nucleic Acids Res.* 15, 1311-1326.
- Felgner, P. L., Holm, M., and Chan, H. (1989). Cationic Liposome Mediated Transfection. *Proc. West. Pharmacol. Soc.* 32, 115-121.
- Felgner, P. L. a., and Ringold, G. M. (1989). Cationic Liposome-Mediated Transfection. *Nature* 337, 387-388.
- Graham, F. L., Smiley, J., Russell, W. C., and Nairn, R. (1977). Characteristics of a Human Cell Line Transformed by DNA from Human Adenovirus Type 5. *J. Gen. Virol.* 36, 59-74.
- Jainchill, J. L., Aaronson, S. A., and Todaro, G. J. (1969). Murine Sarcoma and Leukemia Viruses: Assay Using Clonal Lines of Contact-Inhibited Mouse Cells. *J. Virol.* 4, 549-553.
- Kao, F. T., and Puck, T. T. (1968). Genetics of Somatic Mammalian Cells, VII. Induction and Isolation of Nutritional Mutants in Chinese Hamster Cells. *Proc. Natl. Acad. Sci. USA* 60, 1275-1281.
- Kit, S., Dubbs, D. R., DeTorres, R. A., and Melnick, J. L. (1965). Enhanced Thymidine Kinase Activity Following Infection of Green Monkey Kidney Cells by Simian Adenoviruses, Simian Papovavirus SV40, and an Adenovirus-SV40 "Hybrid". *Virology* 27, 453-457.
- O'Gorman, S., Fox, D. T., and Wahl, G. M. (1991). Recombinase-Mediated Gene Activation and Site-Specific Integration in Mammalian Cells. *Science* 251, 1351-1355.
- Sauer, B. (1994). Site-Specific Recombination: Developments and Applications. *Curr. Opin. Biotechnol.* 5, 521-527.
- Shigekawa, K., and Dower, W. J. (1988). Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. *BioTechniques* 6, 742-751.
- Talavera, A., and Basilico, C. (1977). Temperature Sensitive Mutants of BHK Cells Affected in Cell Cycle Progression. *J. Cell. Physiol.* 92, 425-436.
- Weiss, A., Wiskocil, R. L., and Stobo, J. D. (1984). The Role of T3 Surface Molecules in the Activation of Human T Cells: A Two-stimulus Requirement for IL2 Production Reflects Events Occurring at a Pre-translational Level. *J. Immunol.* 133, 123-128.
- Wigler, M., Silverstein, S., Lee, L.-S., Pellicer, A., Cheng, Y.-C., and Axel, R. (1977). Transfer of Purified Herpes Virus Thymidine Kinase Gene to Cultured Mouse Cells. *Cell* 11, 223-232.

---

©1999-2003 Invitrogen Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.